


# OVAL 8 FINGER SPLINTS

## FEATURES

- Each splint fits two sizes. A size 6 Oval-8 splint fits both a size 6 and a size 6-1/2 depending on which end of the splint is fit on the finger first.
- To immobilize fingers for maximum protection
- Choose the size that fits the IP joint and the next size larger to fit on top.
- Provide effective control at DIP to treat Mallet Finger injuries without limiting function.
- Stabilize the fingers and prevent unwanted posture to facilitate handwriting skills



ITEM #	SIZE
82-P10080-2500	#2 REFILLS 5 PACK
82-P10080-3500	#3 REFILLS 5 PACK
82-P10080-4500	#4 REFILLS 5 PACK
82-P10080-5500	#5 REFILLS 5 PACK
82-P10080-6500	#6 REFILLS 5 PACK
82-P10080-7500	#7 REFILLS 5 PACK
82-P10080-8500	#8 REFILLS 5 PACK
82-P10080-9500	#9 REFILLS 5 PACK
82-P10080-1050	#10 REFILLS 5 PACK
82-P10080-1150	#11 REFILLS 5 PACK
82-P10080-1250	#12 REFILLS 5 PACK
82-P10080-1350	#13 REFILLS 5 PACK
82-P10080-1450	#14 REFILLS 5 PACK
82-P10080-1550	#15 REFILLS 5 PACK