

# TOE LOOP

## FEATURES

- The unique cushioning and grip of our foam-lined material makes these straps ideal for repositioning your toes.
- Loop and wrap design secures two or more toes together to stabilize fracture and derotate toes that turn under.
- Use Toe Loops to cushion the top of your toe and help straighten flexible hammertoes and claw toes.
- The layer of thin, dense foam offers superior pressure distribution without addubg bulk.
- Easy to apply and so light and comfortable, Toe Loops can be worn in any type of footwear.



ITEM #	SIZE
82-P50010-2500	TOE LOOP NARROW 25 PACK
82-P50010-5000	TOE LOOP NARROW 5 PACK
82-P50020-2500	TOE LOOP WIDE 25 PACK
82-P50020-5000	TOE LOOP WIDE 5 PACK